

# RTD™



WOUND DRESSING

4 x 4 x 1/4 in / 10 X 10 x 0.6 cm

17.8x7.6x22.2cm



WOUND DRESSINGS  
4 x 4 x 1/4 in / 10 X 10 x 0.6 cm  
RTD 76040414

RTD™  
WOUND DRESSING

# RTD™

WOUND DRESSING

Size: 4" x 4" x 1/4"  
(10 x 10 x 0.6 cm)

Quantity: (10) Dressing

REF RTD76040414

HCPCS Code: A6209

**STERILE** **Quemada Estéril, no lo utilice si está abierto o dañado antes de la unidad**

INDICATIONS:  
For the treatment of any clinical condition, RTD™  
Wound Dressings are not intended to replace  
absorbent therapy or appropriate course of  
treatment for infection. When dressing on  
wound or burn, do not allow contact with  
any potential irritants which may be present.  
Do not use on open areas. For a complete list of contraindications, see package insert.

Rx Only



FOR MODERATE TO HEAVY EXUDATE  
PARA USO EN HERIDAS CON EXUDADO  
MODERADO A ALTO



Un solo uso

1/4" DRESSING

Manufactured for:



LOT : XX/XXXXX-X  
EXPIRY: YYYY-MM



30856227004045

Keneric Healthcare, LLC, 105 Wilcox Creek Boulevard, Suite 200, Allen, TX 75015, USA - Tel: +1 833 672 2033  
KENERIC HEALTHCARE S.A., Boulevard General Wavre 85, 1020 Brussels, Belgium

10

WOUND DRESSINGS  
4 x 4 x 1/4 in / 10 X 10 x 0.6 cm  
RTD 76040414

10

KENERICHEALTHCARE.COM

B76040414CP-1

10

WOUND DRESSINGS  
4 x 4 x 1/4 in / 10 X 10 x 0.6 cm  
RTD 76040414

RTD™  
WOUND DRESSING

LOT : XX/XXXXX-X  
EXPIRY: YYYY-MM



# RTD™

## WOUND DRESSING

Size: 4" x 4" x 1/4"  
(10 x 10 x 0.6 cm)

Quantity: (1) Dressing

REF RTD71040414

HCP Code: A6209

RTD™ Wound Dressings are comprised of: Absorbent Polyether/Polyurethane; Methylene Blue (up to 0.25 mg/g); Gentian Violet (up to 0.25 mg/g); Silver Sodium Hydrogen Zirconium Phosphate (up to 7mg/g).  
El Apósito para Heridas de RTD™ se compone de espuma de poliuretano/ poliéter, Azul de Metileno (0.25 mg/g), Violeta de Genciana (0.25 mg/g), y Fosfato de Plata, Sodio y Zirconio (7 mg/g).

### PRECAUTIONS

In the event of any clinical infection, RTD™ Wound Dressing is not intended to replace systemic therapy or appropriate course of treatment for infection. When choosing a wound dressing, clinicians should consider any potential adverse effects for any patient. Stop using straight away if an allergic reaction is suspected.

STERILE R

Gamma Estéril; no lo utilice si está abierto o dañado embalaje de la unidad

Rx Only



No echo con Látex



Packaged in China



Consulte las instrucciones



FOR MODERATE TO HEAVY EXUDATE  
PARA USO EN HERIDADS CON EXUDADO  
MODERADO A ALTO



Un solo uso

1/4" DRESSING

Manufactured for:



LOT : XX/XXXXX-X  
EXPIRY: YYYY-MM

 Keneric Healthcare, LLC - 825 Walters Creek Boulevard, Suite 250, Allen, TX 75013, USA - Tel +1 855 872 2013

 Obelis S.A., Boulevard General Wahis 53, 1030 Brussels, Belgium



[www.KenericHealthcare.com](http://www.KenericHealthcare.com)

P71040414CP-1

